


ARK5 Antibody (R35382)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35382-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 9891 |
| Applications | Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000 |
| Limitations | This ARK5 antibody is available for research use only. |



Western blot testing of kidney lysate with ARK5 antibody at 0.3ug/ml. Predicted molecular weight: ~75kDa.

Description

Additional name(s) for this target protein: AMPK-related protein kinase 5; NUAK1

Application Notes

Optimal dilution of the ARK5 antibody should be determined by the researcher.

Immunogen

Amino acids QNRPRPQYLKRYRNR were used as the immunogen for this ARK5 antibody.

Storage

Aliquot and store the ARK5 antibody at -20oC.