


## ARL6 Antibody (R34587)

Catalog No.	Formulation	Size
R34587-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	379
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ARL6 antibody is available for research use only.



Western blot testing of human lung lysate with ARL6 antibody at 0.5ug/ml. Predicted molecular weight: ~22kDa.

## Description

Additional name(s) for this target protein: ADP-ribosylation factor-like 4D/6; ARL4D

## Application Notes

Optimal dilution of the ARL6 antibody should be determined by the researcher.

## Immunogen

Amino acids KQDQPGALSAAEVEKR were used as the immunogen for this ARL6 antibody.

## Storage

Aliquot and store the ARL6 antibody at -20oC.