


Aromatase Antibody / CYP19A1 (R34169)

Catalog No.	Formulation	Size
R34169-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1588
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Aromatase antibody is available for research use only.



Western blot testing of human placenta lysate with Aromatase antibody at 0.1ug/ml.
Predicted molecular weight ~58kDa.

Description

Additional name(s) for this target protein: Cytochrome P450, family 19, subfamily A, polypeptide 1; CYP19A1

Application Notes

Optimal dilution of the Aromatase antibody should be determined by the researcher.

Immunogen

Amino acids HDLSLHPDETKN were used as the immunogen for this Aromatase antibody.

Storage

Aliquot and store the Aromatase antibody at -20oC.