


ARPC4 Antibody (R33961)

Catalog No.	Formulation	Size
R33961-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10093
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ARPC4 antibody is available for research use only.



Western blot testing of human lymph node lysate with ARPC4 antibody at 0.03ug/ml.
Predicted molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: Actin related protein 2/3 complex, subunit 4

Application Notes

Optimal dilution of the ARPC4 antibody should be determined by the researcher.

1. This ARPC4 antibody will detect isoform a of the protein.

Immunogen

Amino acids ERHNKPEVEVR were used as the immunogen for this ARPC4 antibody.

Storage

Aliquot and store the ARPC4 antibody at -20oC.