


Arrestin 3 Antibody / ARR3 (R36442)

Catalog No.	Formulation	Size
R36442-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	407
Localization	Cytoplasmic, membranous
Applications	IHC (FFPE) : 2-4ug/ml
Limitations	This Arrestin 3 antibody is available for research use only.



IHC testing of FFPE human retina with Arrestin 3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining. Staining of the outer segments of rod receptor cells is seen.

Description

Additional name(s) for this target protein: Retinal (x-arrestin); ARR3

Application Notes

Optimal dilution of the Arrestin 3 antibody should be determined by the researcher.

Immunogen

Amino acids RKGEEESQKAVE were used as the immunogen for this Arrestin 3 antibody.

Storage

Aliquot and store the Arrestin 3 antibody at -20oC.