


ARSA Antibody (R36351)

Catalog No.	Formulation	Size
R36351-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	410
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ARSA antibody is available for research use only.



Western blot testing of mouse testis lysate with ARSA antibody at 0.5ug/ml. Expected molecular weight: 53-63 kDa.

Description

Additional name(s) for this target protein: Arylsulfatase A

Application Notes

Optimal dilution of the ARSA antibody should be determined by the researcher.

Immunogen

Amino acids YDLSKDPGENYN were used as the immunogen for this ARSA antibody.

Storage

Aliquot and store the ARSA antibody at -20oC.