


## ARSA Antibody (R36351)

Catalog No.	Formulation	Size
R36351-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	410
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ARSA antibody is available for research use only.



Western blot testing of mouse testis lysate with ARSA antibody at 0.5ug/ml. Expected molecular weight: 53-63 kDa.

## Description

Additional name(s) for this target protein: Arylsulfatase A

## Application Notes

Optimal dilution of the ARSA antibody should be determined by the researcher.

## Immunogen

Amino acids YDLSKDPGENYN were used as the immunogen for this ARSA antibody.

## **Storage**

Aliquot and store the ARSA antibody at -20oC.