


Aryl Hydrocarbon Receptor Antibody (R34997)

Catalog No.	Formulation	Size
R34997-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	196
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Aryl Hydrocarbon Receptor antibody is available for research use only.



IHC testing of FFPE human kidney with Aryl Hydrocarbon Receptor antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: AH receptor; AHR

Visit our [AHR Antibody](#) page to discover additional antibodies against this Environmental Sensor Protein, including reagents for studies of xenobiotic metabolism, immune regulation, and transcriptional signaling.

Application Notes

Optimal dilution of the Aryl Hydrocarbon Receptor antibody should be determined by the researcher.

Immunogen

Amino acids PENQKHGLNPQSA were used as the immunogen for this Aryl Hydrocarbon Receptor antibody.

Storage

Aliquot and store the Aryl Hydrocarbon Receptor antibody at -20oC.