


Arylsulfatase C Antibody / Steroid sulfatase / Steryl-sulfatase (R35191)

Catalog No.	Formulation	Size
R35191-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	412
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Arylsulfatase C antibody is available for research use only.



Western blot testing of human placenta lysate with Arylsulfatase C antibody at 0.01ug/ml. Predicted molecular weight: ~66 kDa.

Description

Additional name(s) for this target protein: Steryl-sulfate sulfohydrolase, STS

Application Notes

Optimal dilution of the Arylsulfatase C antibody should be determined by the researcher.

Immunogen

Amino acids QAGQKIDEPTSN were used as the immunogen for this Arylsulfatase C antibody.

Storage

Aliquot and store the Arylsulfatase C antibody at -20°C.