


Arylsulfatase D Antibody (R33384)

Catalog No.	Formulation	Size
R33384-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	414
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Arylsulfatase D antibody is available for research use only.



Western blot testing of A459 lysate with Arylsulfatase D antibody with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~65/42kDa (isoforms alpha/beta).

Description

Additional name(s) for this target protein: ARSD

Application Notes

Optimal dilution of the Arylsulfatase D antibody should be determined by the researcher.

Immunogen

Amino acids DASNGYRALQWNA were used as the immunogen for this Arylsulfatase D antibody.

Storage

Aliquot and store the Arylsulfatase D antibody at -20°C.