


Ascl3 Antibody (R34231)

Catalog No.	Formulation	Size
R34231-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56676
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Ascl3 antibody is available for research use only.



Western blot testing of mouse liver lysate with Ascl3 antibody at 0.1 ug/ml. Predicted molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: Achaete-scute complex homolog 3

Application Notes

Optimal dilution of the Ascl3 antibody should be determined by the researcher.

1. Not suitable for use in IHC testing.

Immunogen

Amino acids HLPEDYLEKRLSK were used as the immunogen for this Ascl3 antibody.

Storage

Aliquot and store the Ascl3 antibody at -20oC.