

## Ascl3 Antibody (R34231)

Catalog No.	Formulation	Size
R34231-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56676
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Ascl3 antibody is available for research use only.



Western blot testing of mouse liver lysate with Ascl3 antibody at 0.1 ug/ml. Predicted molecular weight: ~20kDa.

## Description

Additional name(s) for this target protein: Achaete-scute complex homolog 3

## Application Notes

Optimal dilution of the Ascl3 antibody should be determined by the researcher.

1. Not suitable for use in IHC testing.

## Immunogen

Amino acids HLPEDYLEKRLSK were used as the immunogen for this Ascl3 antibody.

## Storage

Aliquot and store the Ascl3 antibody at -20oC.