


# ASF1A Antibody (R33901)

Catalog No.	Formulation	Size
R33901-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25842
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ASF1A antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: CCG1-interacting factor A, silencing function 1A

### **Application Notes**

Optimal dilution of the ASF1A antibody should be determined by the researcher.

### **Immunogen**

Amino acids ENSLNVMLESH were used as the immunogen for this ASF1A antibody. **Storage** Aliquot and store the ASF1A antibody at -20oC.