


## ASIC1 Antibody (R34824)

Catalog No.	Formulation	Size
R34824-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	41
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ASIC1 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with ASIC1 antibody at 0.3ug/ml.  
Predicted molecular weight: ~65kDa.

## Description

Additional name(s) for this target protein: Acid-sensing (proton-gated) ion channel 1

## Application Notes

Optimal dilution of the ASIC1 antibody should be determined by the researcher.

## Immunogen

Amino acids NILPHHPARGT were used as the immunogen for this ASIC1 antibody.

## **Storage**

Aliquot and store the ASIC1 antibody at -20oC.