


ASNSD1 Antibody (R35080)

Catalog No.	Formulation	Size
R35080-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54529
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ASNSD1 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human ASNSD1-FLAG probed with ASNSD1 antibody [1ug/ml] in Lane A and anti-FLAG (1/10000) in Lane C. Mock-transfected HEK293 probed with ASNSD1 antibody [1ug/ml] in Lane B. Predicted molecular weight: ~72kDa.

Description

Additional name(s) for this target protein: Asparagine synthetase domain containing 1

Application Notes

Optimal dilution of the ASNSD1 antibody should be determined by the researcher.

Immunogen

Amino acids PNKHVSVSPDR were used as the immunogen for this ASNSD1 antibody.

Storage

Aliquot and store the ASNSD1 antibody at -20oC.