


Aspa Antibody (R36217)

Catalog No.	Formulation	Size
R36217-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11484
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Aspa antibody is available for research use only.



Western blot testing of mouse kidney lysate with Aspa antibody at 0.01ug/ml. Predicted molecular weight: ~35kDa.

Description

Additional name(s) for this target protein: Aspartoacylase; Acy-2

Application Notes

Optimal dilution of the Aspa antibody should be determined by the researcher.

Immunogen

Amino acids TLSAKSIRSTLH were used as the immunogen for this Aspa antibody.

Storage

Aliquot and store the Aspa antibody at -20oC.