


ATF4 Antibody (R33763)

Catalog No.	Formulation	Size
R33763-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	468
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ATF4 antibody is available for research use only.



Western blot testing of K562 lysate with ATF4 antibody at 0.3ug/ml. Predicted molecular weight ~39 kDa.

Description

Additional name(s) for this target protein: Activating transcription factor 4

Application Notes

Optimal dilution of the ATF4 antibody should be determined by the researcher.

Immunogen

Amino acids EEVRKARGKKRVP were used as the immunogen for this ATF4 antibody.

Storage

Aliquot and store the ATF4 antibody at -20oC.