


ATGL Antibody (R34182)

Catalog No.	Formulation	Size
R34182-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57104
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ATGL antibody is available for research use only.



Western blot testing of human heart lysate with ATGL antibody at 0.03ug/ml. Predicted molecular weight ~55kDa.

Description

Additional name(s) for this target protein: Adipose triglyceride lipase, Desnutrin; PNPLA2

Application Notes

Optimal dilution of the ATGL antibody should be determined by the researcher.

Immunogen

Amino acids HFNSKDELIQAN were used as the immunogen for this ATGL antibody.

Storage

Aliquot and store the ATGL antibody at -20oC.