


AURKA Antibody (R35377)

Catalog No.	Formulation	Size
R35377-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6790
Applications	Western Blot : 0.1-0.3ug/ml Immunoprecipitation : suitable under native conditions ELISA (peptide) LOD : 1:128000
Limitations	This AURKA antibody is available for research use only.



Western blot testing of Jurkat lysate with AURKA antibody at 0.1ug/ml. Predicted molecular weight ~45 kDa.

Description

Additional name(s) for this target protein: Aurora kinase A

Application Notes

Optimal dilution of the AURKA antibody should be determined by the researcher.

Immunogen

Amino acids QNKESASKQS were used as the immunogen for this AURKA antibody.

Storage

Aliquot and store the AURKA antibody at -20oC.