


AURKC Antibody (R35129)

Catalog No.	Formulation	Size
R35129-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6795
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This AURKC antibody is available for research use only.



Western blot testing of human testis lysate with AURKC antibody at 1.5ug/ml. Predicted molecular weight: ~36kDa.

Description

Additional name(s) for this target protein: Aurora kinase C, AIE2, AIK3, ARK3, AurC

Application Notes

Optimal dilution of the AURKC antibody should be determined by the researcher.

1. This AURKC antibody is specific for isoform 1.

Immunogen

Amino acids PRAVVQLGKAQP were used as the immunogen for this AURKC antibody.

Storage

Aliquot and store the AURKC antibody at -20oC.