


Aurora B Antibody (R36378)

Catalog No.	Formulation	Size
R36378-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9212
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Aurora B antibody is available for research use only.



Western blot testing of thymus [A], lymph node [B], tonsil [C], spleen [D] and testis [E] lysates with Aurora B antibody at 1ug/ml. Predicted molecular weight: 39-45 kDa

Description

Additional name(s) for this target protein: AURKB; aurora kinase B antibody; AIK2

Application Notes

Optimal dilution of the Aurora B antibody should be determined by the researcher.

Immunogen

Amino acids YKELQKSCTFDEQ were used as the immunogen for this Aurora B antibody.

Storage

Aliquot and store the Aurora B antibody at -20oC.