


## Aurora B Antibody (R36378)

Catalog No.	Formulation	Size
R36378-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9212
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Aurora B antibody is available for research use only.



Western blot testing of thymus [A], lymph node [B], tonsil [C], spleen [D] and testis [E] lysates with Aurora B antibody at 1ug/ml. Predicted molecular weight: 39-45 kDa

## Description

Additional name(s) for this target protein: AURKB; aurora kinase B antibody; AIK2

## Application Notes

Optimal dilution of the Aurora B antibody should be determined by the researcher.

## Immunogen

Amino acids YKELQKSCTFDEQ were used as the immunogen for this Aurora B antibody.

## **Storage**

Aliquot and store the Aurora B antibody at -20oC.