

Autotaxin Antibody (R34975)

Catalog No.	Formulation	Size
R34975-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5168
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Autotaxin antibody is available for research use only.



Western blot testing of human placenta lysate with Autotaxin antibody at 0.3ug/ml.
Expected molecular weight: 99-125 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: ATX, ATX-X, ectonucleotide pyrophosphatase/phosphodiesterase 2; ENPP2

Application Notes

Optimal dilution of the Autotaxin antibody should be determined by the researcher.

Immunogen

Amino acids PDHLTSCVRPDVRVS were used as the immunogen for this Autotaxin antibody.

Storage

Aliquot and store the Autotaxin antibody at -20oC.