


## Autotaxin Antibody (R34975)

Catalog No.	Formulation	Size
R34975-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	5168
<b>Applications</b>	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This Autotaxin antibody is available for research use only.



Western blot testing of human placenta lysate with Autotaxin antibody at 0.3ug/ml.  
Expected molecular weight: 99-125 kDa depending on glycosylation level.

### Description

Additional name(s) for this target protein: ATX, ATX-X, ectonucleotide pyrophosphatase/phosphodiesterase 2; ENPP2

### Application Notes

Optimal dilution of the Autotaxin antibody should be determined by the researcher.

### Immunogen

Amino acids PDHLTSCVRPDVRVS were used as the immunogen for this Autotaxin antibody.

## **Storage**

Aliquot and store the Autotaxin antibody at -20oC.