


AVPR1A Antibody (R33248)

Catalog No.	Formulation	Size
R33248-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	552
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This AVPR1A antibody is available for research use only.



Western blot testing of NIH3T3 lysate with AVPR1A antibody at 0.03ug/ml. Predicted molecular weight ~47/85kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Arginine vasopressin receptor 1A

Application Notes

Optimal dilution of the AVPR1A antibody should be determined by the researcher.

Immunogen

Amino acids CHPLKTLQQPARRSR were used as the immunogen for this AVPR1A antibody.

Storage

Aliquot and store the AVPR1A antibody at -20oC.