


## AVPR1B Antibody (R34558)

Catalog No.	Formulation	Size
R34558-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26361
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AVPR1B antibody is available for research use only.



Western blot testing of human amygdala lysate with AVPR1B antibody, [cat # R35132](#) (0.5ug/ml, lane A) and rat brain lysate with AVPR1B antibody, cat # R34558 (0.05ug/ml, lane B). Predicted molecular weight: ~47/170kDa (unmodified/modified).

## Description

Additional name(s) for this target protein: Arginine vasopressin receptor 1B

## Application Notes

Optimal dilution of the AVPR1B antibody should be determined by the researcher.

## Immunogen

Amino acids KPAGSLKDLEQVD were used as the immunogen for this AVPR1B antibody.

## **Storage**

Aliquot and store the AVPR1B antibody at -20oC.