


## Axotrophin Antibody (R35788)

Catalog No.	Formulation	Size
R35788-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64844
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Axotrophin antibody is available for research use only.



Western blot testing of mouse skeletal muscle lysate with Axotrophin antibody at 2ug/ml.  
Predicted molecular weight: ~76kDa.

## Description

Additional name(s) for this target protein: AXO, AXOT, MARCH7, MARCH-VII

## Application Notes

Optimal dilution of the Axotrophin antibody should be determined by the researcher.

## **Immunogen**

Amino acids RTLQAHMEDLETSED were used as the immunogen for this Axotrophin antibody.

## **Storage**

Aliquot and store the Axotrophin antibody at -20oC.