


B3GNT2 Antibody (R33827)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33827-100UG | 0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|----------------------|---|
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 10678 |
| Applications | Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This B3GNT2 antibody is available for research use only. |



Description

Additional name(s) for this target protein: Beta-1,3-N-acetylglucosaminyltransferase 2

Application Notes

Optimal dilution of the B3GNT2 antibody should be determined by the researcher.

Immunogen

Amino acids EKHKGFRTFDIE were used as the immunogen for this B3GNT2 antibody. **Storage** Aliquot and store the B3GNT2 antibody at -20oC.