


B7-H4 Antibody (R35936)

Catalog No.	Formulation	Size
R35936-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79679
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This B7-H4 antibody is available for research use only.



Western blot testing of human pancreas lysate with B7-H4 antibody at 0.5ug/ml.
Predicted molecular weight ~31kDa.

Description

Additional name(s) for this target protein: B7 family member, H4; VTCN1; V-set domain containing T cell activation inhibitor 1

Application Notes

Optimal dilution of the B7-H4 antibody should be determined by the researcher.

Immunogen

Amino acids SKGKGNNANLEYK were used as the immunogen for this B7-H4 antibody.

Storage

Aliquot and store the B7-H4 antibody at -20oC.