


## BAALC Antibody (R34596)


Catalog No.	Formulation	Size
R34596-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79870
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This BAALC antibody is available for research use only.



Western blot testing of mouse brain lysate with BAALC antibody at 0.03ug/ml. Predicted molecular weight: ~16kDa.



Western blot of HEK293 lysate overexpressing BAALC and probed with BAALC antibody (mock transfection in lane 1). Predicted molecular weight: ~16kDa.

## Description

Additional name(s) for this target protein: Brain and acute leukemia, cytoplasmic

## Application Notes

Optimal dilution of the BAALC antibody should be determined by the researcher.

1. This BAALC antibody is specific for isoform 1.

## Immunogen

Amino acids KQNGQLQTTEAKRDAK were used as the immunogen for this BAALC antibody.

## Storage

Aliquot and store the BAALC antibody at -20oC.