


# **BACH1 Antibody (R34163)**

Catalog No.	Formulation	Size
R34163-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	571
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This BACH1 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: BTB and CNC homology 1

## **Application Notes**

Optimal dilution of the BACH1 antibody should be determined by the researcher.

# **Immunogen**

Amino acids HASVQPNERSENEC were used as the immunogen for this BACH1 antibody.

# **Storage** Aliquot and store the BACH1 antibody at -20oC.