


BAG2 Antibody [clone AEOE-2] (RQ4814)

Catalog No.	Formulation	Size
RQ4814	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AEOE-2
Purity	Affinity purified
UniProt	O95816
Applications	Western Blot : 1:500-1:2000
Limitations	This BAG2 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with BAG2 antibody. Predicted molecular weight ~22 kDa.

Description

BAG family molecular chaperone regulator 2 inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.

Application Notes

Optimal dilution of the BAG2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human BAG2 was used as the immunogen for the BAG2 antibody.

Storage

Store the BAG2 antibody at -20oC.