


## BAG5 Antibody (R33579)

Catalog No.	Formulation	Size
R33579-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9529
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This BAG5 antibody is available for research use only.



Western blot testing of Jurkat lysate with BAG5 antibody at 1ug/ml. Predicted molecular weight: ~56/51kDa (isoforms a/b).

## Description

Additional name(s) for this target protein: BCL2-associated athanogene 5

## Application Notes

Optimal dilution of the BAG5 antibody should be determined by the researcher.

1. This BAG5 antibody is specific for isoforms a and b of the protein.

## **Immunogen**

Amino acids DPQGEEKCKAARKQ were used as the immunogen for this BAG5 antibody.

## **Storage**

Aliquot and store the BAG5 antibody at -20oC.