


BAIAP2 Antibody [isoform 3] (R35979)


Catalog No.	Formulation	Size
R35979-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10458
Localization	Cytoplasmic
Applications	Western Blot : 1-3ug/ml (Transfected) IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This BAIAP2 antibody is available for research use only.



Western blot testing of HEK293 lysate over expressing human BAIAP2-FLAG probed with BAIAP2 antibody (1 ug/ml) in Lane A and with anti-FLAG (1/3000) in lane C. Mock-transfected HEK293 probed with anti-BAIAP2 (1ug/ml) in Lane B. Predicted molecular weight: ~57 kDa.



IHC testing of FFPE human brain (cortex) tissue with BAIAP2 antibody at 3.75ug/ml.
Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: BAI1-associated protein 2

Application Notes

Optimal dilution of the BAIAP2 antibody should be determined by the researcher.

1. This BAIAP2 antibody is specific to isoform 3.

Immunogen

Amino acids SMSSADVEVARF were used as the immunogen for this BAIAP2 antibody.

Storage

Aliquot and store the BAIAP2 antibody at -20oC.