

BAK1 Antibody (R33172)

Catalog No.	Formulation	Size
R33172-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	578
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BAK1 antibody is available for research use only.



Western blot testing of A431 lysate with BAK1 antibody at 0.5ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: BCL2-antagonist/killer 1

Application Notes

Optimal dilution of the BAK1 antibody should be determined by the researcher.

Immunogen

Amino acids ASGQGPGPPRQE were used as the immunogen for this BAK1 antibody.

Storage

Aliquot and store the BAK1 antibody at -20oC.