


BAT3 Antibody / BAG6 (R35001)

Catalog No.	Formulation	Size
R35001-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7917
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This BAT3 antibody is available for research use only.



IHC testing of FFPE human spleen (nucleoli staining) with BAT3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: HLA-B associated transcript 3, BAG6, BAG-6; BAT3

Application Notes

Optimal dilution of the BAT3 antibody should be determined by the researcher.

Immunogen

Amino acids PETNAPNHPSAEY were used as the immunogen for this BAT3 antibody.

Storage

Aliquot and store the BAT3 antibody at -20oC.