


Bax Antibody / BCL2-associated X (R32761)


Catalog No.	Formulation	Size
R32761	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q07812
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml Flow Cytometry : 1-3ug/million cells Immunofluorescence : 5ug/ml
Limitations	This Bax antibody is available for research use only.




IHC staining of FFPE human lung cancer tissue with Bax antibody, HRP-labeled secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE mouse kidney tissue with Bax antibody, HRP-labeled secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE rat kidney tissue with Bax antibody, HRP-labeled secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Immunofluorescent staining of FFPE human A549 cells with Bax antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) human ThP-1, 2) human MCF7, 3) human A549, 4) human HeLa and 5) rat C6 cell lysate with Bax antibody at 0.5ug/ml. Predicted molecular weight: 21-24 kDa.



Flow cytometry testing of fixed and permeabilized human ThP-1 cells with Bax antibody at 1ug/10⁶ cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= Bax antibody.

Description

Bax Antibody specifically binds Bax, also known as BCL2-associated X protein, a pro-apoptotic member of the BCL2 protein family that plays a critical role in the intrinsic pathway of apoptosis. It promotes programmed cell death by permeabilizing the mitochondrial outer membrane, leading to the release of cytochrome c and subsequent activation of caspases. Bax is tightly regulated at the transcriptional and post-translational levels, and its activity is modulated through interactions with other BCL2 family proteins.

Bax exists primarily in the cytosol in healthy cells but undergoes conformational changes and translocates to the mitochondria upon apoptotic stimuli. It forms homo-oligomers that are essential for pore formation and mitochondrial

membrane disruption. This makes Bax a key regulator of cellular homeostasis, development, and response to stress or damage.

Multiple isoforms of Bax have been reported, resulting from alternative splicing. These isoforms may differ in their pro-apoptotic potency, subcellular localization, and expression profiles, making it important to use a Bax antibody that recognizes the relevant forms in specific experimental contexts.

A high-performance Bax antibody is essential for studies involving apoptosis, cancer, neurodegeneration, and immune response. The Bax antibody is widely used in applications such as western blot, immunohistochemistry, immunofluorescence, and flow cytometry. A well-validated Bax antibody can reliably detect expression levels and cellular distribution of Bax across different models and conditions.

NSJ Bioreagents offers a range of Bax antibody products designed for sensitivity, specificity, and versatility across research applications. Each Bax antibody is quality tested to ensure reproducible results in your apoptosis-related experiments.

Learn more about Bax biology and explore our highly published mouse monoclonal Bax antibody on our [Bax Antibody](#) page.

Application Notes

Optimal dilution of the Bax antibody should be determined by the researcher.

Immunogen

Amino acids 17-48 (EQIMKTGALLLQGFIQDRAGRMGGEAPELALD) from the human protein were used as the immunogen for the Bax antibody.

Storage

After reconstitution, the Bax antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.