


## BBS4 Antibody (R34785)

Catalog No.	Formulation	Size
R34785-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	585
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BBS4 antibody is available for research use only.



Western blot testing of HeLa lysate with BBS4 antibody at 1ug/ml. Predicted molecular weight: ~58kDa.

## Description

Additional name(s) for this target protein: Bardet-Biedl syndrome 4; BBS4

## Application Notes

Optimal dilution of the BBS4 antibody should be determined by the researcher.

## Immunogen

Amino acids NEAAKLNQKDWEISH were used as the immunogen for this BBS4 antibody.

## **Storage**

Aliquot and store the BBS4 antibody at -20oC.