


BCAR3 Antibody (R35653)

Catalog No.	Formulation	Size
R35653-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8412
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BCAR3 antibody is available for research use only.



Description

Additional name(s) for this target protein: Breast cancer estrogen resistance 3

Application Notes

Optimal dilution of the BCAR3 antibody should be determined by the researcher.

Immunogen

Amino acids RKLEPPPVKQAEL were used as the immunogen for this BCAR3 antibody.

Storage

Aliquot and store the BCAR3 antibody at -20oC.