


BCAR3 Antibody (R35653)

Catalog No.	Formulation	Size
R35653-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8412
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BCAR3 antibody is available for research use only.



Western blot testing of mouse kidney epithelial cells lysate with BCAR3 antibody 0.5ug/ml. Predicted molecular weight: ~93kDa, observed here at ~110kDa.

Description

Additional name(s) for this target protein: Breast cancer estrogen resistance 3

Application Notes

Optimal dilution of the BCAR3 antibody should be determined by the researcher.

Immunogen

Amino acids RKLEPPPVKQAEI were used as the immunogen for this BCAR3 antibody.

Storage

Aliquot and store the BCAR3 antibody at -20°C.