


BCL2 Antibody (R33135)

Catalog No.	Formulation	Size
R33135-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	596
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This BCL2 antibody is available for research use only.



Western blot testing of human lymph lysate with BCL2 antibody at 0.1ug/ml. Expected molecular weight: ~26kDa (alpha) and ~22kDa (beta).

Description

Additional name(s) for this target protein: B-cell CLL/lymphoma 2

Application Notes

Optimal dilution of the BCL2 antibody should be determined by the researcher.

Immunogen

Amino acids AGRTGYDNREIVMKYC were used as the immunogen for this BCL2 antibody.

Storage

Aliquot and store the BCL2 antibody at -20oC.