


BCL7A Antibody (R34752)

Catalog No.	Formulation	Size
R34752-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	605
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This BCL7A antibody is available for research use only.



Western blot testing of human lymph node lysate with BCL7A antibody at 1ug/ml.
Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: B-cell CLL/lymphoma 7A

Application Notes

Optimal dilution of the BCL7A antibody should be determined by the researcher.

Immunogen

Amino acids MKLEASQQNSEEM were used as the immunogen for this BCL7A antibody.

Storage

Aliquot and store the BCL7A antibody at -20°C.