

BCOR Antibody (R34741)

Catalog No.	Formulation	Size
R34741-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54880
Applications	ChIP : suitable ELISA (peptide) LOD : 1:2000
Limitations	This BCOR antibody is available for research use only.



Description

Additional name(s) for this target protein: BCL-6 interacting corepressor

Application Notes

Optimal dilution of the BCOR antibody should be determined by the researcher.

Immunogen

Amino acids LSATPLYGNVHSW were used as the immunogen for this BCOR antibody.

Storage

Aliquot and store the BCOR antibody at -20oC.