

BDKRB1 Antibody (R34695)

Catalog No.	Formulation	Size
R34695-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	623
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This BDKRB1 antibody is available for research use only.



Western blot testing of K562 lysate with BDKRB1 antibody at 2ug/ml. Predicted molecular weight: ~40kDa, but can be observed at ~35kDa.

Description

Additional name(s) for this target protein: Bradykinin receptor B1

Application Notes

Optimal dilution of the BDKRB1 antibody should be determined by the researcher.

Immunogen

Amino acids KVVELYKQCTPK were used as the immunogen for this BDKRB1 antibody.

Storage

Aliquot and store the BDKRB1 antibody at -20oC.