


## BDNF Antibody (R34016)

Catalog No.	Formulation	Size
R34016-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	627
Applications	Immunofluorescence : 10ug/ml
Limitations	This BDNF antibody is available for research use only.



Immunofluorescent staining of fixed and permeabilized human U-251 cells with BDNF antibody (green) at 10ug/ml and DAPI nuclear stain.

## Description

Additional name(s) for this target protein: Brain-derived neurotrophic factor

## Application Notes

Optimal dilution of the BDNF antibody should be determined by the researcher.

## Immunogen

Amino acids ETKCNPMGYTKE were used as the immunogen for this BDNF antibody.

## Storage

Aliquot and store the BDNF antibody at -20oC.