


BDNF Antibody (R34016)

Catalog No.	Formulation	Size
R34016-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	627
Applications	Immunofluorescence : 10ug/ml
Limitations	This BDNF antibody is available for research use only.



Immunofluorescent staining of fixed and permeabilized human U-251 cells with BDNF antibody (green) at 10ug/ml and DAPI nuclear stain.

Description

Additional name(s) for this target protein: Brain-derived neurotrophic factor

Application Notes

Optimal dilution of the BDNF antibody should be determined by the researcher.

Immunogen

Amino acids ETKCNPMGYTKE were used as the immunogen for this BDNF antibody.

Storage

Aliquot and store the BDNF antibody at -20oC.