

## Beta-2-Microglobulin Antibody (R35447)

Catalog No.	Formulation	Size
R35447-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	567
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Beta-2-Microglobulin antibody is available for research use only.



Western blot testing of human duodenum lysate with Beta-2-Microglobulin antibody at 0.3ug/ml. Predicted molecular weight ~14 kDa.

## Description

Additional name(s) for this target protein: B2M

## Application Notes

Optimal dilution of the Beta-2-Microglobulin antibody should be determined by the researcher.

## **Immunogen**

Amino acids QRTPKIQVYSRHP were used as the immunogen for this Beta-2-Microglobulin antibody.

## **Storage**

Aliquot and store the Beta-2-Microglobulin antibody at -20oC.