


BID Antibody (R36322)

Catalog No.	Formulation	Size
R36322-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	637
Localization	Cytoplasmic
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This BID antibody is available for research use only.



Western blot testing of A431 lysate with BID antibody at 0.1ug/ml. Predicted molecular weight ~26/22 kDa (isoforms 1/2)

Description

Additional name(s) for this target protein: BH3 interacting domain death agonist

Application Notes

Optimal dilution of the BID antibody should be determined by the researcher.

Immunogen

Amino acids TYVRSLARNGMD were used as the immunogen for this BID antibody.

Storage

Aliquot and store the BID antibody at -20oC.