


## Biglycan Antibody (R33614)

Catalog No.	Formulation	Size
R33614-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Rat
Predicted Reactivity	Dog, Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	633
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Biglycan antibody is available for research use only.



Western blot testing of rat skin lysate with Biglycan antibody at 0.1ug/ml. Predicted molecular weight: ~42 kDa.

## Description

Additional name(s) for this target protein: Bone/cartilage proteoglycan-I; BGN

## Application Notes

Optimal dilution of the Biglycan antibody should be determined by the researcher.

## Immunogen

Amino acids DRLAIQFGNYKK were used as the immunogen for this Biglycan antibody.

## Storage

Aliquot and store the Biglycan antibody at -20oC.