


BIRC3 Antibody (R36028)

Catalog No.	Formulation	Size
R36028-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	330
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BIRC3 antibody is available for research use only.



Western blot testing of HeLa lysate with BIRC3 antibody at 0.1ug/ml. Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: Baculoviral IAP repeat-containing 3

Application Notes

Optimal dilution of the BIRC3 antibody should be determined by the researcher.

Immunogen

Amino acids SQPTFPSSVT were used as the immunogen for this BIRC3 antibody.

Storage

Aliquot and store the BIRC3 antibody at -20oC.